



LUND UNIVERSITY

Multi-elemental Analysis of Aerosol Samples with PIXE in Combination with Complementary Nuclear Techniques

Bohgard, Mats; Akselsson, Roland

1982

[Link to publication](#)

Citation for published version (APA):

Bohgard, M., & Akselsson, R. (1982). *Multi-elemental Analysis of Aerosol Samples with PIXE in Combination with Complementary Nuclear Techniques*. Abstract from 10th Conference of Gesellschaft für Aerosolforschung, Bologna, Italy.

Total number of authors:

2

General rights

Unless other specific re-use rights are stated the following general rights apply:

Copyright and moral rights for the publications made accessible in the public portal are retained by the authors and/or other copyright owners and it is a condition of accessing publications that users recognise and abide by the legal requirements associated with these rights.

- Users may download and print one copy of any publication from the public portal for the purpose of private study or research.
- You may not further distribute the material or use it for any profit-making activity or commercial gain
- You may freely distribute the URL identifying the publication in the public portal

Read more about Creative commons licenses: <https://creativecommons.org/licenses/>

Take down policy

If you believe that this document breaches copyright please contact us providing details, and we will remove access to the work immediately and investigate your claim.

LUND UNIVERSITY

PO Box 117
221 00 Lund
+46 46-222 00 00

Gaet, Bologna, Italy, 1982

MULTIELEMENTAL ANALYSIS OF AEROSOL SAMPLES WITH PIXE IN COMBINATION
WITH COMPLEMENTARY NUCLEAR TECHNIQUES 37

Mats Bohgard and K Roland Akselsson 27

Department of Nuclear Physics, Lund Institute of Technology, Sölvegatan 14, S-223 62 Lund, Sweden

During the last decade a new analytical method, Particle Induced X-ray Emission Analysis (PIXE), has been developed. The method has been extensively applied to aerosol studies. The capability of the method for reliable multielemental analysis with low detection limits (10^{-9} - 10^{-12} g) for several elements of interest, has been attractive. The possibility to analyse small samples ($\sim 10^{-6}$ g) has given opportunities for using small and low-weight equipment for size- and/or time-fractionated sampling.

In PIXE-analysis the samples are bombarded with ions from an accelerator and the emitted characteristic X-rays are detected with an energy-dispersive semiconductor detector. The energy of the characteristic X-rays decreases with decreasing atomic number of the elements. For the lightest elements (Li, Be, B, C, N, O, F, Mg, Al, Si and P) either the X-ray energy is too low for detection or the attenuation of the X-ray intensity in a thin aerosol sample significantly increases the uncertainty of the analysis.

However, complementary nuclear techniques may be used to determine the light elements. Particle Elastic Scattering analysis (PESA) and detection of radiation from particle-induced nuclear reactions are two techniques which can be exploited. Sometimes the complementary techniques may be performed simultaneously with the PIXE-analysis.

For optimum use of PIXE and complementary methods the sampling technique and the choice of sampling substrate as well as the analysis parameters should be considered. The potential of PIXE and complementary methods will be discussed and examples of applications of the combination sampling technique, PIXE analysis and complementary techniques will be given.